

# 10 Reasons to Choose ESR-Auto Plus®

An automated sedimentation rate analyzer.



MORE THAN 2,600 UNITS PLACED



- **>98% Correlation to Modified Westergren Method**  
Provides accurate results in only 30 minutes (15-minute prediction mode). Ensures confidence and saves time and money.
- **ESR Placement Program**  
The automated instrument is available for purchase or through an ESR Placement Program.
- **Internal Printer**  
Automatically prints patient results at the conclusion of the test. Eliminates repeat testing and saves valuable benchtop space.
- **Full Random Access**  
Eliminates sample batching and provides flexibility. Ten-sample capacity allows patient testing at your convenience.
- **Positive Sample Identification**  
Barcode scanner aids in patient sample identification and prevents transcription errors. Stores up to 500 patient results and 100 QC results per level in log file.
- **LIS Compatible**  
An LIS interface prevents transcription errors and aids in patient sample identification.
- **72-Hour Patient Sample Stability**  
ESR-Vacuum Tubes preserve the integrity of the patient sample for up to 72 hours at 2 °C to 10 °C. Ensures specimen stability during transit from off-site blood collection facility.
- **Safe Blood Collection System**  
Safety-coated tubes reduce exposure to potentially hazardous material. All ESR-Vacuum Tubes allow direct vacuum blood collection or transfer from an EDTA (purple-top) tube.
- **Manual Back-Up System**  
The ESR-10 Manual Rack produces results in 30 minutes and uses Streck ESR-Vacuum Tubes and ESR-Chex control.
- **Excellent Service**  
Streck provides superior sales, customer and technical service to help you through the evaluation and ordering process. The ESR-Auto Plus has a proven track record and substantial peer group for ILQC (inter-laboratory quality control) reporting.

**Streck**

800.843.0912  
streck.com